

Nucleus Segmentation Cell zone estimation 3DSpatialTissueJ plugin

<https://github.com/nhuhoa/Spatial3DTissueJ>

Hoa Tran

Workflow for tissue spatial organization analysis




Figure : Workflow: (i) segmentation and cell phenotype detection (ii) spatial organization analysis

3D Nuclei Segmentation




Figure : Steps of nuclei segmentation.




Figure : Diagram of nuclei segmentation [1].


Nucleus Segmentation Plugin

Noted: in windows system, your path will be:
ex: C::\yourdir\inputfolder\

Using '\' instead of '/' as in linux sys


Browser button automatically take care of directory path syntax

- Example macro
- 3 parameters need to defined:
- `dir="C:/Users/SALAB VR/Documents/Hoa/Spatial3DTissueJ-master/
H1536_islet1/";`
- **`small_nucleus_diameter_thrs=15;`**
- **`large_nucleus_diameter_thrs=36;`**
- `open(dir+"C4-dapi.tif");`
- ***/// Using automatic seeds selection mode***
- `run("NUCLEI SEGMENTATION", "save_dir=["+dir+"] dapi=C4-dapi.tif
maximal="+large_nucleus_diameter_thrs+"
minimal="+small_nucleus_diameter_thrs+" seed=29000 automatic");`



Parameter 1: seed threshold

- Seeds: local maximal, pixels with intensity higher than intensity of the neighbourhood pixels.
- Manually define seed threshold value: ex: 29000, 30000 for 16 bits (similar to ex:120 /255 for 8 bits image)
- Automatic define seed threshold: (program will generate a seed threshold value)
 - Get mean value of all pixels in image
 - Seed = mean + sd, ex: sd=500, a pixel is considered as a seed if its intensity is higher than mean value of all pixels.




Parameter 2: nucleus diameters

- Example macro
- `dir="C:/Users/SALAB VR/Documents/Hoa/Spatial3DTissueJ-master/H1536_islet1/";`
- **`small_nucleus_diameter_thrs=15;`**
- **`large_nucleus_diameter_thrs=36;`**
- `open(dir+"C4-dapi.tif");`
- *//// Using automatic seeds selection mode*
- `run("NUCLEI SEGMENTATION", "save_dir=["+dir+"] dapi=C4-dapi.tif maximal="+large_nucleus_diameter_thrs+" minimal="+small_nucleus_diameter_thrs+" seed=29000 automatic");`



Key idea: Fourier transform as the transformation function for seed detection.

Small nuc
diameter
threshold



$$H(x, y, z) = \begin{cases} 1 & \text{if } D(x, y, z) \leq D_0 \\ 0 & \text{otherwise} \end{cases}$$

D_{max}





remove a big
part of noise




$$H(x, y, z) = \begin{cases} 1 & \text{if } D(x, y, z) \geq D_1 \\ 0 & \text{otherwise} \end{cases}$$

D_{min}





highlight edge info
and sharpen



$$H(x, y, z) = \begin{cases} 1 & \text{if } D(x, y, z) > D_0 \& D(x, y, z) < D_1 \\ 0 & \text{if } D(x, y, z) < D_0 \cup D(x, y, z) > D_1 \end{cases}$$

$[D_{min} - D_{max}]$



combine the feature
of low and high pass




Large nuc
diameter
threshold

Figure : Band pass filtering and implementation in our 3D nuclei segmentation application.

Parameter 2: nucleus diameters

In current version, I define the nucleus diameters by measuring nucleus using rule tool.

In the future version, I can improve program by automatic estimate nuc radius but need time to test new algo






- To measure the small and large nucleus radius diameters: first we convert micron to pixels scale:
- ImageJ menu → Image → Properties
- Convert scale from microns to pixel

Parameter 2: nucleus diameters


Using ruler tool to measure nucleus diameter

- Using ruler tool
- Look at several nucleus: a big one, a small one, and check nucleus diameters.
- You can measure one time and use this setting parameter for other islet images in case you use same microscopy setting (ex: same microns to pixels) for all captured images
 - *Ex: in plugin, I used: small_nucleus_diameter_thrs=15;*
 - *large_nucleus_diameter_thrs=36;*



Parameter 2: nucleus diameters

- Look at the demo_only.tif image, if you see the nucleus are detected as your expectation —> good output. Otherwise, you can tune diameters thresholds by increasing or decreasing diameters thresholds.



Workflow for tissue spatial organization analysis




Figure : Workflow: (i) segmentation and cell phenotype detection (ii) spatial organization analysis

Cell zone estimation

- Cell zone estimation: from nucleus
 - red area, extend into space to obtain cell zone area.
- In general: extend as much as possible —> not good
- Improved here: define the maximum of extension: ex: N pixels from nucleus border





Figure : Cell zone estimation using watershed algorithms.


Cell zone estimation

- Example macro:
- ```
dir="C:/Users/SALAB VR/Documents/
Hoa/Spatial3DTissueJ-master/
H1536_islet1/";
```
- ```
max_radius_extension_wat=5; //change
this threshold if your cell zone area in the
image is smaller, depend on image
resolution, in general 3 to 5 is a good one.
```
- ```
open(dir+"C4-dapi.tif");
```
- ```
open(dir+"dapi-seg.tif");
```
- ```
run("CELL ZONE ESTIMATION",
"save_dir=["+dir+"] nuclei=dapi-seg.tif
radius_max="+max_radius_extension_wat
+" save");
```



# Parameter: radius of extension

- Using ruler tool
- Look at several nucleus: and check the maximum extension from nucleus border to the cell periphery
- You can measure one time and use this setting parameter for other islet images in case you use same microscopy setting (ex: same microns to pixels) for all captured images
  - *Ex: in plugin, I used: max\_radius\_extension\_wat=5;*



# Acknowledgements

- [1] "A novel toolbox to investigate tissue spatial organization applied to the study of the islets of Langerhans"  
Hoa Tran et al.
- [2] J. Ollion, J. Cochenne, F. Loll, C. Escudé, T. Boudier. (2013) TANGO: A Generic Tool for High-throughput 3D Image Analysis for Studying Nuclear Organization. *Bioinformatics* 2013 Jul 15;29(14):1840-1.
- [3] Spatial3DTissueJ plugins would like to thank P. Andrey, J.-F. Gilles and the developers of the following plugins :
  - ImageScience
  - LocalThickness
  - ConvexHull3D
  - 3D Object Counter
  - Droplet Counter